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Read carefully prior to use

WARNING


This product contains a chemical (Ethylene Oxide) known to the State of California to cause cancer, birth defects or other reproductive harm.

CAUTION: FEDERAL (U.S.A.) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

 **INTEGRA™**
Neuro Sciences
CAMINO®

MICRO VENTRICULAR BOLT PRESSURE-TEMPERATURE MONITORING KIT MODEL 110-4HMT

Figure 1

SYSTEM DESCRIPTION

The Camino® Micro Ventricular Pressure-Temperature Monitoring Kit consists of a sterile transducer-tipped pressure monitoring catheter with thermistor and accessory items which provide a means of access to the cerebral ventricles for CSF sampling and drainage, fluid injection and continuously monitoring intracranial pressure and temperature. Since its method of measuring pressure is unique, please read this section carefully.

Unlike ordinary pressure monitoring systems, the Camino Catheter has a miniature transducer and thermistor at the distal tip. This unique design eliminates the need for a "fluid-filled system" to carry pressure waves to an external transducer. Thus, the problems inherent in such fluid-filled systems are eliminated. The result is a precise pressure measurement and an artifact free, high fidelity waveform trace.

The Micro Ventricular Pressure-Temperature Monitoring Kit contains the following accessory items for use with the Camino Catheter (Figure 1):

- Camino Ventricular Access Device with Ventricular Catheter, Bolt, Female Luer Lock and Stylet,
- 0.209 inch (5.31 mm) diameter No. 4 Drill Bit with Safety Stop,
- Additional Luer Lock Port to connect drainage system,
- Hex Wrench,
- Zero Adjustment Tool,

- Thermistor Connector Cover.

Note: A hand drill and various standard surgical instruments and supplies are required to place the ICP/ICT monitoring system.

A set of instruments and supplies is available from Integra NeuroSciences in our Cranial Access Kits.

INDICATIONS

The Camino Micro Ventricular Bolt Pressure-Temperature Monitoring Kit is indicated for use by qualified neurosurgeons for measurement of intracranial pressure and temperature in the ventricles and for cerebrospinal fluid drainage. The Camino Micro Ventricular Bolt Pressure-Temperature Monitoring Kit is intended to be used with an external drainage system as indicated by individual manufacturers.

CONTRAINDICATIONS

This device is not intended for any use other than that indicated. This device is contraindicated for use in the MRI field.

RISK AND COMPLICATIONS

- Maintenance of sterility during placement and subsequent handling of the Camino Micro Ventricular Bolt Pressure-Temperature Monitoring Kit is essential. Sterile technique should be used at all times when inserting, adjusting, and securing the Camino Catheter.
- If monitoring is continued for more than 5 days, placement of a new system under sterile conditions is recommended.
- Infection, subcutaneous leakage of CSF, neurological sequelae, and blockage by intraventricular debris (including bloody and/or highly proteinic CSF) have occurred during the use of ventricular catheters.
- Placement of the tip opening within the reach of choroid plexus has resulted in blockage of ventricular catheters.
- The ventriculostomy must be carried out by a qualified neurosurgeon using standard surgical procedures and skill. Hemorrhage from the dura or cortical surface at the bolt insertion site may occur. Patients should be tested for normal blood clotting function prior to bolt placement. Appropriate measures to avoid infections and complications are the sole responsibility of the neurosurgeon in charge.
- This device utilizes PVC tubing which is known to contain DEHP in a portion of the drainage path which on occasion is used for sampling and fluid injection. The risk of exposure to the male fetus, male neonate, and peri-pubescent male is not considered significant for the following reasons: Sampled and drained CSF is not intended to be reintroduced into the patient. Injected fluids such as anti-coagulants and saline solution are only used as needed to eliminate or reduce blockage of the ventricular catheter and thus are not continuous in use.

PRECAUTIONS

- It is imperative that the ventricular catheter not be handled with bare fingers or come into contact with linty surfaces. Silicone elastomers are very electrostatic and therefore susceptible to contamination by airborne or surface particles. The presence of these contaminants could cause adverse tissue reaction. Rubber-shod clamps or washed, gloved hands are the best means of handling implantable silicone devices.
- Extreme bending and/or kinks can impair the performance of the Fiber Optic Pressure Transducer. Exercise caution when handling the catheter.
- The catheter is designed for **SINGLE USE ONLY. DO NOT RESTERILIZE OR REUSE.** Camino 110-4HMT is supplied in a single use package and is guaranteed to be sterile

and non-pyrogenic unless opened or damaged. The product is intended for use as an intracranial pressure-temperature monitor and is not to be reused. Any attempts to resterilize or reuse the product will damage the catheter and impair its ability to function as intended.

- The in-dwelling portion of the ventricular catheter assembly is fabricated of barium-impregnated silicone elastomer. Its OD is 3.7mm and ID is 2.2 mm. Its length is adjustable from 6-8 cm. Designed for directing fluid from the ventricles through a series of drainage holes, it can be inserted into the ventricular cavity with the supplied stainless steel stylet.
- Use aseptic technique throughout procedures.
- Maintain the insertion site with regular meticulous redressing using aseptic technique.
- Do not attach anything to transducer air vent. Vent must remain open for proper operation (Figure 1).

CAUTION

- Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- Connect the thermistor connector only to Host Monitors marked "Patient connection electrically isolated" or "BF" or "CF" or marked with the international symbols:

or

- When connecting the thermistor connector to a Host Monitor, refer to the Host Monitor's operations manual for complete instructions.
- Verify proper operation of the combined Host Monitor-110-4HMT System before clinical use.
- The combined leakage currents of devices interconnected with the 110-4HMT can lead to a potentially hazardous condition. When combined with other devices/systems, ensure that the combined system leakage current does not exceed 0.1 mA.

INSTRUCTIONS FOR USE

INSERTION METHOD FOR THE CAMINO MICRO VENTRICULAR BOLT

The recommended frontal placement is 3-4 cm off the midline, just anterior to the coronal suture. After the site has been chosen, the area is shaved and prepped in a sterile fashion, usually with a betadine solution. The shaved and prepared area is then draped. The area of

the incision is infiltrated subcutaneously with 1% Xylocaine. An approximately three centimeter linear incision is made and carried to the bone. A self-retaining retractor is then inserted to provide good bone exposure and hemostasis of the skin edges.

Adjust the safety stop on the drill bit to the estimated skull thickness and secure firmly with the hex wrench.

Secure the drill bit to a twist drill and in a standard fashion drill a hole through the outer and inner tables of the skull, taking care to minimize any potential for parenchymal injury. Penetrate the dura under direct vision with a #11 blade, securing hemostasis as necessary.

Using the stylet, insert the ventricular catheter into the ventricle. When the CSF is obtained, hold the catheter securely, remove the stylet, slide the bolt down and screw in, using bone wax to ensure a tight seal. Do not over tighten, as stripping of the threads may cause loss of seal.

Continue to hold the catheter securely, and turn the compression cap clockwise to lock the catheter in place. Slide the strain relief down and attach to the compression cap. Cap the catheter with Luer cap to prevent CSF loss.

CAMINO PRESSURE-TEMPERATURE MONITORING CATHETER PREPARATION PRIOR TO INSERTION INTO VENTRICULAR CATHETER

The Camino Catheter can be used with the Camino 420/V420 series or with MPM/SPM series pressure monitors. For Camino monitor set-up and use, refer to the appropriate Camino monitor Directions For Use.


Remove the Camino Catheter from its sterile package and firmly attach the connector to the pre-amp connector (Figure 2). If the monitor display does not read zero after a short system self-check delay, use the tool from the catheter kit to turn the zero adjustment on the bottom side of the transducer connector until the monitor display reads zero (Figure 3).

Remove cover from the thermistor connector of the Camino Catheter. To interface with host monitor cardiac output module, connect cardiac output monitor cable connector to Camino Thermistor Connector. Verify that a temperature is displayed on the Host Monitor.

INSERTION OF CAMINO CATHETER

Remove the Luer cap from the ventricular catheter, insert the Camino Transducer-tipped Catheter and secure Luer lock. Holding the ventricular catheter straight will facilitate passage.

Prepare an external ventricular drainage system according to its manufacturer's directions and

attach to the side port of the Y-connector. The user should close the drainage system for approximately five minutes before recording the intracranial pressure (ICP) during monitoring. This will allow the system to adjust for under-reading the ICP values caused by air, drainage bag leveled below the foramen of Monro, gravity effect on the tubing that may be hanging off the bed, and the potential rise in the ICP that occurs when the system is closed off with the stopcocks closest to the patient's head. The system would be treated as a "trending" value if the stopcock is not closed to record the value and the user should recognize that the value could be lower than the closed value.

When monitoring is to be discontinued, detach the strain relief from the compression cap. Loosen compression cap and remove ventricular catheter prior to the removal of the bolt from the skull.

IT IS RECOMMENDED THAT THE CATHETER BE DISCONNECTED FROM THE PREAMP CONNECTOR WHEN THE PATIENT IS MOVED. THIS WILL NOT AFFECT CALIBRATION.

CONTINUOUS PRESSURE AND TEMPERATURE MONITORING

Since the Camino Catheter has a miniaturized transducer at the distal tip, it requires no fluid-filled system. Thus, the need for an external transducer, pressure dome, and pressure tubing is eliminated. As a result, temperature and pressure may be monitored continuously without flushing or recalibration.

PRODUCT INFORMATION DISCLOSURE
INTEGRA NEUROSCIENCES HAS EXERCISED


REASONABLE CARE IN THE SELECTION OF MATERIALS AND THE MANUFACTURE OF THESE PRODUCTS. INTEGRA NEUROSCIENCES EXCLUDES ALL WARRANTIES, WHETHER EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. INTEGRA NEUROSCIENCES SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM USE OF THIS PRODUCT. INTEGRA NEUROSCIENCES NEITHER ASSUMES NOR AUTHORIZES ANY PERSON TO ASSUME FOR IT ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THESE PRODUCTS.

RETURNED GOODS POLICY


Products must be returned in unopened packages with manufacturer's seals intact to be accepted for replacement or credit unless returned due to a complaint of product defect. Determination of a product defect will be made by Integra NeuroSciences. Products will not be accepted for replacement if they have been in the possession of the customer for more than 90 days.

**RECOMMENDED
INSERTION
PROCEDURE
SUMMARY**


After drilling hole in skull and penetrating the dura,


1. Pass ventricular catheter assembly into ventricle.


2. Remove stylet.


3. Slide bolt down and screw in.


4. Secure compression cap.


5. Slide strain relief sheath down and attach to compression cap.

6. *Zero the transducer*, then insert transducer-tipped catheter into ventricular catheter.

7. Secure luer lock.

8. Connect drainage system to y-site.

INTERNATIONAL MARKINGS

CATHETER EXTENSION CABLE CONNECTIONS

Symbol also signifies:

TYPE BF EQUIPMENT: Protected against electric shock, having an F-type isolated (floating) applied part.

Sterilized using ethylene oxide.

Date of expiration.

Do not reuse, single use device.

Manufacturer.

Do not resterilize.

Not for Use in an MRI field.

TYPE CF EQUIPMENT: Provided with a degree of protection against electric shock higher than Type BF equipment, having an F-type isolated (floating) applied part; intended for direct cardiac application.

Refer to Directions for Use.

Lot Number.

Rx ONLY

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Authorized representative in the European Community.

Contains or presence of DEHP.

SPECIFICATIONS

Ventricular Catheter:	
Outside Diameter	3.7 mm
Inside Diameter	2.2 mm
Length	6-8 cm (adjustable)
Transducer size	4F (1.35mm Diameter)
Transducer type	Fiber optic
Frequency response (system)	
Model MPM-1/SPM-1	100Hz (-3dB)
Model V420/M420/420XP	120Hz (-3dB)
Model 420	33Hz (-3dB)
Measurement range (system)	-10 to 125 mmHg
Zero drift (system)	
First 24 hours (maximum)	0 ± 2 mmHg
5 days (typical)	less than ± 1 mmHg per day
Reference pressure	Atmosphere
Overpressure	-700 to 1250 mmHg
Temperature coefficient	Max of 3 mmHg over temperature range of 22°C to 38°C (70°F-100°F)
Linearity and hysteresis (system)	
Pressure range:	
-10 to 50 mmHg	± 2 mmHg or better
51 to 125 mmHg	± 6% of reading or better
Temperature Measurement	
Temperature Measurement Resolution	0.1°C
Temperature Measurement Range	30°C - 40°C
Temperature Measurement Accuracy	± 0.3°C (30°C - 40°C)

Camino and Integra NeuroSciences are registered trademarks at Integra LifeSciences Corporation in the United States and/or other countries. The Integra wave logo is a trademark of Integra LifeSciences Corporation © 2010 Integra LifeSciences Corporation. All Rights Reserved.

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Integra NeuroSciences Ltd.,